

Methionine aminopeptidase 2

Catalog No: #AP80628



Package Size: #AP80628-1 50ug #AP80628-2 100ug #AP80628-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Methionine aminopeptidase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MAP 2, MetAP 2, Initiation factor 2-associated 67 kDa glycoprotein, Peptidase M 2, p67eIF2, p67, Mnpep, P67eif2
Accession No.	Uniprot ID:O08663
Uniprot	O08663
GeneID	56307;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Metap2

Note: This product is for in vitro research use only