

Protein arginine N-methyltransferase 1

Catalog No: #AP80629



Package Size: #AP80629-1 50ug #AP80629-2 100ug #AP80629-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein arginine N-methyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Histone-arginine N-methyltransferase PRMT1,Hrmt112,Mrmt1
Accession No.	Uniprot ID:Q9JIF0
Uniprot	Q9JIF0
GeneID	15469;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Prmt1

Note: This product is for in vitro research use only