

Apolipoprotein B

Catalog No: #AP80631



Package Size: #AP80631-1 50ug #AP80631-2 100ug #AP80631-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Apolipoprotein B |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Chicken |
| Accession No. | Uniprot ID:P11682 |
| Uniprot | P11682 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: APOB

Note: This product is for in vitro research use only