

Cyclic AMP-dependent transcription factor ATF-1

Catalog No: #AP80632



Package Size: #AP80632-1 50ug #AP80632-2 100ug #AP80632-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclic AMP-dependent transcription factor ATF-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	cAMP-dependent transcription factor ATF-1,Activating transcription factor 1
Accession No.	Uniprot ID:Q08DA8
Uniprot	Q08DA8
GeneID	506967;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ATF1

Note: This product is for in vitro research use only