

RAC-alpha serine/threonine-protein kinase

Catalog No: #AP80636



Package Size: #AP80636-1 50ug #AP80636-2 100ug #AP80636-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RAC-alpha serine/threonine-protein kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AKT1 kinase,Protein kinase B,PKB,Protein kinase B alpha,PKB alpha,Proto-oncogene c-Akt,RAC-PK-alpha,Thymoma viral proto-oncogene,Akt,Rac
Accession No.	Uniprot ID:P31750
Uniprot	P31750
GeneID	11651;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Akt1

Note: This product is for in vitro research use only