

AKT-interacting protein

Catalog No: #AP80638



Package Size: #AP80638-1 50ug #AP80638-2 100ug #AP80638-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AKT-interacting protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Fused toes protein homolog,FTS,RCJMB04_17k13
Accession No.	Uniprot ID:Q5ZJJ5
Uniprot	Q5ZJJ5
GeneID	415720;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AKTIP

Note: This product is for in vitro research use only