

Aldehyde dehydrogenase, mitochondrial

Catalog No: #AP80639



Package Size: #AP80639-1 50ug #AP80639-2 100ug #AP80639-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aldehyde dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AHD-M1,ALDH class 2,ALDH-E2,ALDHI,Ahd-1,Ahd1
Accession No.	Uniprot ID:P47738
Uniprot	P47738
GeneID	11669;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aldh2

Note: This product is for in vitro research use only