

Fructose-bisphosphate aldolase B

Catalog No: #AP80640



Package Size: #AP80640-1 50ug #AP80640-2 100ug #AP80640-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fructose-bisphosphate aldolase B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Liver-type aldolase,ALDB
Accession No.	Uniprot ID:P05062
Uniprot	P05062
GeneID	229;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ALDOB

Note: This product is for in vitro research use only