

Fructose-bisphosphate aldolase B

Catalog No: #AP80641



Package Size: #AP80641-1 50ug #AP80641-2 100ug #AP80641-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fructose-bisphosphate aldolase B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Aldolase 2,Liver-type aldolase,Aldo2
Accession No.	Uniprot ID:Q91Y97
Uniprot	Q91Y97
GeneID	230163;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aldob

Note: This product is for in vitro research use only