

Angiotensin-converting enzyme 3

Catalog No: #AP80645



Package Size: #AP80645-1 50ug #AP80645-2 100ug #AP80645-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiotensin-converting enzyme 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Angiotensin-converting enzyme 3, ANGPT3, ANGPT3-like protein 3, ANGPT3L3, UNQ153/PRO179
Accession No.	Uniprot ID:Q9Y5C1
Uniprot	Q9Y5C1
GeneID	27329;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ANGPTL3

Note: This product is for in vitro research use only