

Apoptosis-associated speck-like protein containing a CARD

Catalog No: #AP80647

Package Size: #AP80647-1 50ug #AP80647-2 100ug #AP80647-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Apoptosis-associated speck-like protein containing a CARD
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	mASC,PYD and CARD domain-containing protein,Asc
Accession No.	Uniprot ID:Q9EPB4
Uniprot	Q9EPB4
GeneID	66824;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pycard

Note: This product is for in vitro research use only