

## Angiotensin-converting enzyme 2

Catalog No: #AP80649



Package Size: #AP80649-1 50ug #AP80649-2 100ug #AP80649-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Angiotensin-converting enzyme 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Angiotensin-converting enzyme 2, Heparin-binding EGF-like protein, HEGP
Accession No.	Uniprot ID:Q6TMA8
Uniprot	Q6TMA8
GeneID	362850;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Angptl4

Note: This product is for in vitro research use only