

Disintegrin and metalloproteinase domain-containing protein 8

Catalog No: #AP80650

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80650-1 50ug #AP80650-2 100ug #AP80650-3 1mg

Description

Product Name	Disintegrin and metalloproteinase domain-containing protein 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ADAM 8,Cell surface antigen MS2,Macrophage cysteine-rich glycoprotein,CD156a,Ms2
Accession No.	Uniprot ID:Q05910
Uniprot	Q05910
GeneID	11501;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Adam8

Note: This product is for in vitro research use only