

# Aquaporin-4

Catalog No: #AP80653



Package Size: #AP80653-1 50ug #AP80653-2 100ug #AP80653-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Aquaporin-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AQP-4,Mercurial-insensitive water channel,MIWC,WCH4
Accession No.	Uniprot ID:P55088
Uniprot	P55088
GeneID	11829;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Aqp4

Note: This product is for in vitro research use only