

Annexin A3

Catalog No: #AP80655



Package Size: #AP80655-1 50ug #AP80655-2 100ug #AP80655-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Annexin A3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	35-alpha calcimedlin,Annexin III,Annexin-3,Lipocortin III,Placental anticoagulant protein III,PAP-III,Anx3
Accession No.	Uniprot ID:O35639
Uniprot	O35639
GeneID	11745;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Anxa3

Note: This product is for in vitro research use only