

Arginase-1

Catalog No: #AP80656



Package Size: #AP80656-1 50ug #AP80656-2 100ug #AP80656-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Arginase-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Liver-type arginase, Type I arginase
Accession No.	Uniprot ID:Q95JC8
Uniprot	Q95JC8
GeneID	397115;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ARG1

Note: This product is for in vitro research use only