

Annexin A4

Catalog No: #AP80657



Package Size: #AP80657-1 50ug #AP80657-2 100ug #AP80657-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Annexin A4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Annexin IV,Annexin-4,Anx4
Accession No.	Uniprot ID:P97429
Uniprot	P97429
GeneID	11746;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Anxa4

Note: This product is for in vitro research use only