

Transferrin Monoclonal Antibody(6D6)

Catalog No: #44124

Package Size: #44124-1 50ul #44124-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Transferrin Monoclonal Antibody(6D6)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	6D6
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu
Specificity	Transferrin Mouse Monoclonal AntibodyB£B`6D6B£B© detects endogenous levels of Transferrin
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	Transferrin
Accession No.	Swiss-Prot:P02787Gene ID:7018
Uniprot	P02787
GeneID	7018;
Calculated MW	77kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded human-liver using antibody diluted at 1:50.

Note: This product is for in vitro research use only