

Amphiregulin

Catalog No: #AP80669



Package Size: #AP80669-1 50ug #AP80669-2 100ug #AP80669-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Amphiregulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AR,Colorectum cell-derived growth factor,CRDGF,AREGB,SDGF
Accession No.	Uniprot ID:P15514
Uniprot	P15514
GeneID	374;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AREG

Note: This product is for in vitro research use only