

Apolipoprotein C-III

Catalog No: #AP80675



Package Size: #AP80675-1 50ug #AP80675-2 100ug #AP80675-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Apolipoprotein C-III
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Apo-CIII,Apolipoprotein C3
Accession No.	Uniprot ID:P33622
Uniprot	P33622
GeneID	11814;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Apoc3

Note: This product is for in vitro research use only