

Angiogenin-2

Catalog No: #AP80679



Package Size: #AP80679-1 50ug #AP80679-2 100ug #AP80679-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Angiogenin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Angiogenin-related protein,Angrp
Accession No.	Uniprot ID:Q64438
Uniprot	Q64438
GeneID	11731;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ang2

Note: This product is for in vitro research use only