

Beta-arrestin-1

Catalog No: #AP80680



Package Size: #AP80680-1 50ug #AP80680-2 100ug #AP80680-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Beta-arrestin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Arrestin beta-1,ARR1
Accession No.	Uniprot ID:P49407
Uniprot	P49407
GeneID	408;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ARRB1

Note: This product is for in vitro research use only