

Activity-regulated cytoskeleton-associated protein

Catalog No: #AP80682



Package Size: #AP80682-1 50ug #AP80682-2 100ug #AP80682-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Activity-regulated cytoskeleton-associated protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ARC/ARG3.1,mArc,Activity-regulated gene 3.1 protein homolog,Arg3.1
Accession No.	Uniprot ID:Q9WV31
Uniprot	Q9WV31
GeneID	11838;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Arc

Note: This product is for in vitro research use only