

Antithrombin-III

Catalog No: #AP80692



Package Size: #AP80692-1 50ug #AP80692-2 100ug #AP80692-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Antithrombin-III
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ATIII,Serpin C1,AT3
Accession No.	Uniprot ID:P41361
Uniprot	P41361
GeneID	540261;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SERPINC1

Note: This product is for in vitro research use only