

Beta-2-glycoprotein 1

Catalog No: #AP80698



Package Size: #AP80698-1 50ug #AP80698-2 100ug #AP80698-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Beta-2-glycoprotein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Apolipoprotein H,Apo-H,Beta-2-glycoprotein I,B2GPI,Beta(2)GPI
Accession No.	Uniprot ID:P33703
Uniprot	P33703
GeneID	403945;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:APOH

Note: This product is for in vitro research use only