

## Aprataxin and PNK-like factor

Catalog No: #AP80701



Package Size: #AP80701-1 50ug #AP80701-2 100ug #AP80701-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Aprataxin and PNK-like factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Apurinic-apyrimidinic endonuclease APLF,PNK and APTX-like FHA domain-containing protein,XRCC1-interacting protein 1,C2orf13,PALF,XIP1
Accession No.	Uniprot ID:Q8IW19
Uniprot	Q8IW19
GeneID	200558;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:APLF

Note: This product is for in vitro research use only