

## Isoaspartyl peptidase/L-asparaginase

Catalog No: #AP80703



Package Size: #AP80703-1 50ug #AP80703-2 100ug #AP80703-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Isoaspartyl peptidase/L-asparaginase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Asparaginase-like protein 1, Beta-aspartyl-peptidase, Isoaspartyl dipeptidase, L-asparagine amidohydrolase
Accession No.	Uniprot ID:Q8C0M9
Uniprot	Q8C0M9
GeneID	66514;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Asrg1

Note: This product is for in vitro research use only