

Acetylserotonin O-methyltransferase

Catalog No: #AP80704



Package Size: #AP80704-1 50ug #AP80704-2 100ug #AP80704-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Acetylserotonin O-methyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Hydroxyindole O-methyltransferase,HIOMT
Accession No.	Uniprot ID:P46597
Uniprot	P46597
GeneID	438;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ASMT

Note: This product is for in vitro research use only