

Aprataxin

Catalog No: #AP80709



Package Size: #AP80709-1 50ug #AP80709-2 100ug #AP80709-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Aprataxin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Forkhead-associated domain histidine triad-like protein,FHA-HIT
Accession No.	Uniprot ID:Q8K4H4
Uniprot	Q8K4H4
GeneID	259271;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aptx

Note: This product is for in vitro research use only