

PPAR Delta Mouse Monoclonal Antibody(1D7)

Catalog No: #44132

Package Size: #44132-1 50ul #44132-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PPAR Delta Mouse Monoclonal Antibody(1D7)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1D7
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	PPAR Delta Mouse Monoclonal Antibody(1D7) detects endogenous levels of PPAR Delta
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	PPAR Delta
Other Names	PPAR Delta
Accession No.	Swiss-Prot:Q03181Gene ID:5467
Uniprot	Q03181
GeneID	5467;
Calculated MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using PPAR Delta Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using PPAR Delta Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only