

Amphiregulin

Catalog No: #AP80712



Package Size: #AP80712-1 50ug #AP80712-2 100ug #AP80712-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Amphiregulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	AR,Schwannoma-derived growth factor,SDGF,SdGF
Accession No.	Uniprot ID:P24338
Uniprot	P24338
GeneID	29183;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Areg

Note: This product is for in vitro research use only