

Sarcoplasmic/endoplasmic reticulum calcium ATPase

1



Catalog No: #AP80718

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80718-1 50ug #AP80718-2 100ug #AP80718-3 1mg

Description

Product Name	Sarcoplasmic/endoplasmic reticulum calcium ATPase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	SERCA1, Calcium pump 1, Calcium-transporting ATPase sarcoplasmic reticulum type, fast twitch skeletal muscle isoform, Endoplasmic reticulum class 1/2 Ca(2+) ATPase
Accession No.	Uniprot ID:Q0VCY0
Uniprot	Q0VCY0
GeneID	518117;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ATP2A1

Note: This product is for in vitro research use only