

# A disintegrin and metalloproteinase with thrombospondin motifs 2

Catalog No: #AP80721

Package Size: #AP80721-1 50ug #AP80721-2 100ug #AP80721-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	A disintegrin and metalloproteinase with thrombospondin motifs 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ADAM-TS 2,Procollagen I N-proteinase,PC I-NP,Procollagen I/II amino propeptide-processing enzyme,Procollagen N-endopeptidase,pNPI,PCINP,PCPNI
Accession No.	Uniprot ID:O95450
Uniprot	O95450
GeneID	9509;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ADAMTS2

Note: This product is for in vitro research use only