

Ataxin-10

Catalog No: #AP80722



Package Size: #AP80722-1 50ug #AP80722-2 100ug #AP80722-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ataxin-10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Spinocerebellar ataxia type 10 protein homolog
Accession No.	Uniprot ID:Q2TBW0
Uniprot	Q2TBW0
GeneID	767990;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ATXN10

Note: This product is for in vitro research use only