

## Acrosomal protein SP-10

Catalog No: #AP80724



Package Size: #AP80724-1 50ug #AP80724-2 100ug #AP80724-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Acrosomal protein SP-10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Acrosomal vesicle protein 1,MSA-63
Accession No.	Uniprot ID:P50289
Uniprot	P50289
GeneID	11451;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Acrv1

Note: This product is for in vitro research use only