

Aflatoxin B1 aldehyde reductase member 2

Catalog No: #AP80729



Package Size: #AP80729-1 50ug #AP80729-2 100ug #AP80729-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aflatoxin B1 aldehyde reductase member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	rAFAR2,Succinic semialdehyde reductase,SSA reductase,Afar2,Aiar
Accession No.	Uniprot ID:Q8CG45
Uniprot	Q8CG45
GeneID	171445;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Akr7a2

Note: This product is for in vitro research use only