

Aquaporin-2

Catalog No: #AP80731



Package Size: #AP80731-1 50ug #AP80731-2 100ug #AP80731-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Aquaporin-2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | AQP-2,ADH water channel,Aquaporin-CD,AQP-CD,Collecting duct water channel protein,WCH-CD,Water channel protein for renal collecting duct |
| Accession No. | Uniprot ID:P79099 |
| Uniprot | P79099 |
| GeneID | 539870; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AQP2

Note: This product is for in vitro research use only