

B-cell lymphoma/leukemia 11A

Catalog No: #AP80732



Package Size: #AP80732-1 50ug #AP80732-2 100ug #AP80732-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | B-cell lymphoma/leukemia 11A |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | BCL-11A,B-cell CLL/lymphoma 11A,COUP-TF-interacting protein 1,Ecotropic viral integration site 9 protein,EVI-9,Ctip1,Evi9,Kiaa1809 |
| Accession No. | Uniprot ID:Q9QYE3 |
| Uniprot | Q9QYE3 |
| GeneID | 14025; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Bcl11a

Note: This product is for in vitro research use only