

Baculoviral IAP repeat-containing protein 5

Catalog No: #AP80734



Package Size: #AP80734-1 50ug #AP80734-2 100ug #AP80734-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Baculoviral IAP repeat-containing protein 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Apoptosis inhibitor survivin
Accession No.	Uniprot ID:Q6J1J1
Uniprot	Q6J1J1
GeneID	414925;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: BIRC5

Note: This product is for in vitro research use only