

AT-rich interactive domain-containing protein 4B

Catalog No: #AP80736



Package Size: #AP80736-1 50ug #AP80736-2 100ug #AP80736-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AT-rich interactive domain-containing protein 4B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ARID domain-containing protein 4B, 180 kDa Sin3-associated polypeptide, Sin3-associated polypeptide p180, Histone deacetylase complex subunit SAP180, Retinoblastoma-binding protein 1-like 1, Rbbp111, Sap180
Accession No.	Uniprot ID:Q9JKB5
Uniprot	Q9JKB5
GeneID	84481;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Arid4b

Note: This product is for in vitro research use only