

## JAK1 Mouse Monoclonal Antibody(7G6)

Catalog No: #44135

Package Size: #44135-1 50ul #44135-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	JAK1 Mouse Monoclonal Antibody(7G6)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	7G6
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	JAK1 Mouse Monoclonal Antibody(7G6) detects endogenous levels of JAK1
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	JAK1
Other Names	JAK1
Accession No.	Swiss-Prot:P23458Gene ID:3716
Uniprot	P23458
GeneID	3716;
Calculated MW	130kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using JAK1 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using JAK1 Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only