

Annexin A2

Catalog No: #AP80745



Package Size: #AP80745-1 50ug #AP80745-2 100ug #AP80745-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Annexin A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Annexin II,Annexin-2,Calpactin I heavy chain,Calpactin-1 heavy chain,Chromobindin-8,Lipocortin II,Placental anticoagulant protein IV,PAP-IV,Protein I,p36,Anx2,Cal1h
Accession No.	Uniprot ID:P07356
Uniprot	P07356
GeneID	12306;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Anxa2

Note: This product is for in vitro research use only