

Glutamyl aminopeptidase

Catalog No: #AP80746



Package Size: #AP80746-1 50ug #AP80746-2 100ug #AP80746-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Glutamyl aminopeptidase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | EAP,Aminopeptidase A,AP-A,BP-1/6C3 antigen,CD249 |
| Accession No. | Uniprot ID:P16406 |
| Uniprot | P16406 |
| GeneID | 13809; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Enpep

Note: This product is for in vitro research use only