

Cytosol aminopeptidase

Catalog No: #AP80747



Package Size: #AP80747-1 50ug #AP80747-2 100ug #AP80747-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytosol aminopeptidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Leucine aminopeptidase 3,LAP-3,Leucyl aminopeptidase,Proline aminopeptidase,Prolyl aminopeptidase
Accession No.	Uniprot ID:P28839
Uniprot	P28839
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LAP3

Note: This product is for in vitro research use only