

## Angiotensin-converting enzyme 1

Catalog No: #AP80749

Package Size: #AP80749-1 50ug #AP80749-2 100ug #AP80749-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Angiotensin-converting enzyme 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Angiotensin-converting enzyme 1
Accession No.	Uniprot ID:Q640P2
Uniprot	Q640P2
GeneID	72713;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Angpt1

Note: This product is for in vitro research use only