

Alpha-hemoglobin-stabilizing protein

Catalog No: #AP80751



Package Size: #AP80751-1 50ug #AP80751-2 100ug #AP80751-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha-hemoglobin-stabilizing protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Erythroid-associated factor, ERAF
Accession No.	Uniprot ID:Q865F8
Uniprot	Q865F8
GeneID	338381;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: AHSP

Note: This product is for in vitro research use only