

Interferon-inducible protein AIM2

Catalog No: #AP80752



Package Size: #AP80752-1 50ug #AP80752-2 100ug #AP80752-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interferon-inducible protein AIM2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Interferon-inducible protein 210, Ii-210, Interferon-inducible protein p210, Gm1313, Ii210
Accession No.	Uniprot ID:Q91VJ1
Uniprot	Q91VJ1
GeneID	383619;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Aim2

Note: This product is for in vitro research use only