

A-kinase anchor protein 17A

Catalog No: #AP80753



Package Size: #AP80753-1 50ug #AP80753-2 100ug #AP80753-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	A-kinase anchor protein 17A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AKAP-17A,721P,B-lymphocyte antigen,Protein XE7,Protein kinase A-anchoring protein 17A,PRKA17A,Splicing factor, arginine/serine-rich 17A,CXYorf3,DXYS155E,SFRS17A,XE7
Accession No.	Uniprot ID:Q02040
Uniprot	Q02040
GeneID	8227;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AKAP17A

Note: This product is for in vitro research use only