

RAC-alpha serine/threonine-protein kinase

Catalog No: #AP80754



Package Size: #AP80754-1 50ug #AP80754-2 100ug #AP80754-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	RAC-alpha serine/threonine-protein kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein kinase B,PKB,Protein kinase B alpha,PKB alpha,Proto-oncogene c-Akt,RAC-PK-alpha,PKB,RAC
Accession No.	Uniprot ID:P31749
Uniprot	P31749
GeneID	207;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AKT1

Note: This product is for in vitro research use only