

Proline-rich AKT1 substrate 1

Catalog No: #AP80755



Package Size: #AP80755-1 50ug #AP80755-2 100ug #AP80755-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proline-rich AKT1 substrate 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	40 kDa proline-rich AKT substrate,PRAS40
Accession No.	Uniprot ID:Q96B36
Uniprot	Q96B36
GeneID	84335;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: AKT1S1

Note: This product is for in vitro research use only